



RBD Instruments, Inc.  
2437 NE Twin Knolls Dr.  
Suite 2`  
Bend, Oregon, USA

+1 541-330-0723

<https://rbdinstruments.com/products/vb.html>

## VB FAQ

[Is the length of the IRB-600 Emitter customizable?](#)

Yes. However, there are limitations. Customization requests will be evaluated on a case-by-case basis with respect to feasibility.

[What is the average lifetime of the IRB-600 emitter?](#)

The average lifetime is 5,000 hours. This can be affected by how hot the emitter is run and how close it is to other surfaces (in terms of reflected energy).