



# microCMA Dimensions Form

Name: \_\_\_\_\_

Company: \_\_\_\_\_

Date: \_\_\_\_\_

Comments: \_\_\_\_\_

Signature: \_\_\_\_\_

*By signing this I agree that the dimensions entered on this form meet my specifications.*

Flange to target distance \_\_\_\_\_ mm

+ Z-translator fully compressed length \_\_\_\_\_ 66 mm (50 mm stroke) \_\_\_\_\_ 70 mm (75 mm stroke) \_\_\_\_\_ 74 mm (100 mm stroke)

\_\_\_\_\_ Custom \_\_\_\_\_ mm

- End of optics to target 4 mm

= microCMA length \_\_\_\_\_ mm (standard length 283 mm)

