



microCMA Dimensions Form

Name:	_____
Company:	_____
Date:	_____
Comments:	_____
Signature:	_____

By signing this I agree that the dimensions entered on this form meet my specifications.

Flange to target distance: _____ mm

+ (plus) Z-translator fully compressed length:

	66 mm (50 mm stroke; P/N EQA-BLT27S-2)
	70 mm (75 mm stroke; P/N EQA-BLT27S-3)
	74 mm (100 mm stroke; P/N EQA-BLT27S-4)
	81 mm (150 mm stroke; P/N EQA-BLT27S-6)

- (minus) End of optics to target: **4 mm**

= microCMA length _____ mm (standard length 283 mm)

