



# microCMA Dimensions Form

Name: \_\_\_\_\_

Company: \_\_\_\_\_

Date: \_\_\_\_\_

Comments: \_\_\_\_\_

Signature: \_\_\_\_\_

*By signing this I agree that the dimensions entered on this form meet my specifications.*

Flange to target distance \_\_\_\_\_ mm

+ Z-translator fully compressed length

66 mm  
(50 mm stroke)

70 mm  
(75 mm stroke)

74 mm  
(100 mm stroke)

Custom \_\_\_\_\_ mm

- End of optics to target 4 mm

= microCMA length \_\_\_\_\_ mm (standard length 283 mm)

Vacuum Chamber

