



microCMA Dimensions Form

Name:	_____
Company:	_____
Date:	_____
Comments:	_____
Signature:	_____

By signing this I agree that the dimensions entered on this form meet my specifications.

Flange to target distance: _____ mm

+ (plus) Z-translator fully compressed length:

66 mm	(50 mm stroke; P/N EQA-BLT27S-2)
70 mm	(75 mm stroke; P/N EQA-BLT27S-3)
74 mm	(100 mm stroke; P/N EQA-BLT27S-4)
81 mm	(150 mm stroke; P/N EQA-BLT27S-6)

- (minus) End of optics to target: 4 mm

= microCMA length _____ mm (standard length 283 mm)

