



description

With the exchange of your existing 04-725, 04-727, or 07-763 sample intro arm, RBD can provide you with a magnetic load-lock system at a fraction of the cost of other magnetic load-lock upgrades.

Before



After



benefits

- Better vacuum during sample loading
- System won't vent if the probe is bumped
- No Teflon® seals and spacers to replace
- No expensive intro switch assembly to break
- 2.75" CF Mounting
- Magnetic switch compatible with Perkin-Elmer Auto Valve Controller
- Compatible with standard sample mounts

**PHI is a registered trademark of Physical Electronics
All specifications are subject to change without notice.*